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INSTRUCTIONS FOR USE AND SAFETY

Callus Pro Implants™ (CP)

1. Instructions for use and safety, important information

All products of Callus Implant Systems are recommended exclusively for dentists and oral surgeons who are skilled in dental implantology and oral surgery treatments. In every case, it is the responsibility of the user to decide about the proper utilization of Callus Implant Systems in accordance with the proper medical examination. It is also the responsibility of the user of Callus Implant Systems to decide whether the products are appropriate for a certain patient in a given case. Callus Ltd. disclaims any obligations or responsibility for any direct or indirect damage caused by the implantation or utilization of Callus products among conditions that resulted from decisions which were not appropriately made or which were the results of inexperience of the user. All elements of this surgical and prosthetic set are part of an integrated system, therefore only this set can be utilized safely and successfully – in accordance with instructions and guides included by the supplier – with its original components and devices. In case any part of these sets are used along with any foreign device, it may negatively influence the results, which exempts Callus Implant Solutions Ltd. from any guarantee or other commitments involving legal procedures. Every medical devices and tools within Callus Implant Systems have been developed for special applications. Any improperly utilization may cause these devices to wear off sooner than normal and thus endangers the success of the application.

Drills, devices and all other components of Callus Implant Systems have been developed for different types and diameters. Utilization with inappropriate types or with inappropriate diameters may cause mechanical damage to the components of the system or even tissue damage in the patient. There is a defined surgical set for each implant type. Color marks help with the selection of the proper devices and components. Studying the details of the usage of Callus Implant Systems is the responsibility of the user. In case of any doubt or dubious situations, the user is obliged to contact with Callus Implant Solutions Ltd.

I. IMPLANTS

2. Features and details of Callus Implants

Callus implants are made of high purity titanium designed for use in humans (and Grade 4 quality titanium). The surface of the implants is first sandblasted with large granule bioceramics, and then passivated, namely etched with acid. Abutments and assemblies are made of Grade 5 Ti6Al4V ELI titanium compounds, zirconium oxide ceramics and polymer materials.

3. Indications

Callus Implant System implants can be utilized to provide implantation based/anchored prosthesis for the functional or esthetic rehabilitation of both the lower and the upper jaw. Implants can be used in case of single missing tooth or even to treat partial or even complete edentulous status. Implant abutments can be used to achieve fix prostheses, meaning to fix the bridges, crowns to the implants via cementing or screwing. Alternatively abutments could also enhance removable prostheses/dentures. All indications of various types of Callus implants and head parts will be as instructors for special use as can be found in the guidelines supplied by Callus Implant Systems.

Additional indications

Callus Pro™ Implants can be implanted during surgical exploration intraosseally, under the gum (closing the gum above the implant during the healing time period) or trans gingivally (keeping the gum open with a gum shading during the healing time period). The above procedure is called one- or two-phases surgery. For any further details of implantation see Callus manuals.

4. Contraindications

The most common contraindications of dental implant treatment are the following: recurring mucosal conditions, inflammations of the jaws and gingiva, poor or insufficient oral hygiene, dryness of the mouth (xerostomia), untreated parafunctional habits (grinding, lockjaw), smoking, and still growing jaws (typically under the age of 18).

Dental implant treatment is also contraindicated in case of any alterations or conditions that may complicate surgery, including neurological and psychological diseases, disorders or any other condition that may hinder the patient's practice of regular oral hygiene. Additionally, dental implant treatment is contraindicated in case of any bone alterations due to metabolic disorders (bone metabolic disorders, any disease influencing bone regeneration or micro-circulation, rheumatic disease), untreated hemophilia, or tumors, and 1,5 – 2 years after cancer treatment. Other contraindicating systemic diseases include blood dyscrasia, haemophilia, leukemia, untreated diabetes, hyperthyroidism, AIDS, drug abuse, excessive alcohol consumption or osteoporosis, serious systemic diseases, decreased immunity that increases infection risk, leukocyte dysfunction, endocrine disease, anticoagulation treatment, hemorrhagic diathesis, arteriosclerosis, apoplexy, hypertension, hepatitis, diabetes mellitus. Implantations are also contraindicated in case the patient had a myocardial infarction in the previous 12 months, or in case the patient is immunocompromised or has leukocytosis, or taking steroids regularly, having titanium allergy, or untreated endocrine disorders.

Inappropriate bone quality, radix relicta in the jaw, bone disease or insufficient bone and soft tissue quantity, healing difficulties, oral infections, serious therapy resistant functional disorder, Immunosuppressive treatment, connective tissue diseases/collagenosis, malignoma.

Unfinished bone growth, anatomically insufficient bone status, jaw bone radiation therapy, temporo-mandibular joint disease (TMJ dysfunction), untreated jaw diseases.

Dental implant treatment is contraindicated for patients who treated by bisphosphonate and/ or denosumab containing medicines for 2 years after the medication is stopped. In the time period from 2 to 10 years thorough evaluation should be taken by an Oral Surgeon Specialist evaluating the risk-benefit ratio of the dental implant treatment on the given patient.

Important information:

3.3 and 3.75mm implants in diameter should not be loaded to make up for the single upper first incisors, or any single canine, premolar or molar tooth.

5. Side effects

After surgical implantation, patients should refrain from any physical stress and smoking. Medical instructions and the prescribed therapy must be followed carefully!

Temporary complaints may include pain, facial swelling, speech disorders, gingivitis, hematoma or late onset bleeding. Long-term complaints may include permanent pain caused by the dental implant, permanent paresthesia, dyesthesia, plate peeling, loss of bone tissue from the maxillary/mandibular ridge, local or systemic infections /oroantral or oronasal inflammations/, loss of neighboring teeth, irreversible damage to the neighboring teeth, breaking of the implant, prosthesis, jaw or bone, esthetic problems, local or systemic allergic response. Some cases osseointegration could be failed.

It is the responsibility of the dentist to select patients who are expected to be satisfied psychologically, functionally and esthetically after recovery.

6. Warning

Dental implant surgery and implantation prosthetics are complex dental procedures. It is especially recommended to master all necessary elements of proficiency. An improper surgical procedure may cause implanted loss or bone absorption. The implantation must be preceded by extensive diagnostic procedures. Improper design may result in unsuccessful implantation. Callus implants must only be loaded with instruments included in the Callus Implant Systems and in accordance with guidelines and recommendations supplied therewith.

Getting acquainted with the instructions of usage alone is not enough to start using Callus Implant Systems. Callus implants may exclusively be loaded by oral surgeon specialists or genetal dentists trained to use Callus Systems. Callus Implant Solutions Ltd. is organizing proper in depth training give the user deeper knowledge on this field. Callus Guarantees shall be submitted only if the claimer is able to confirm the completion of a dental implant treatment course organized by Callus Implant Solutions Ltd.

7. Sterilization, packing and storage

Implants are released following gamma sterilization and made for a single use only. Clinically infected im-

plants must not be cleaned and re-sterilized. Callus Implant Solutions Ltd. does not take any liability for the utilization of Callus implants taken from open, re-sterilized Callus implant sets. Used or unsterile implants must not be implanted to bone. The sterile packing protects Callus implants from external contamination. Primary packing must only be opened directly before placing the implant into the bone. Prior to opening the surgeon should check that the packing is not damaged in a way that the implant sterility is at risk. When the Callus Implant is stored properly (at a dry place, protected from direct sunlight at room temperature and in its original packing) implants remain sterile till the use-by date indicated on the external packaging. Implants must not be used after the date of expiry!
All Callus implants can be traceable according to reference number (REF) and batch no.(LOT), which can be found on both the external and the internal packaging's.

8. Handling the implants

Callus *titanium* implants are packed in triple packaging, thus the surface of the implant is protected from mechanical damage or contamination. After removing the inner protective packaging, it is strictly forbidden to touch the implants by hand or any other way. Primarily use the Callus implant surgical screw is defined in the surgical protocol!
If necessary the osseointegration surface of the implant could be touched with sterile surgical device, primarily with tweezer made from pure titanium only!

PREPARING THE IMPLANTS

Each and every Callus implant is packed in a triple-layer sterile packaging. Before use, implant should be stored at room temperature, at a dry place, protected from direct sunlight. Before opening, check the use-by date and the damage if any on the packaging. Packaging should be opened directly before use only. Implant must not be used after the expiration (use-by) date, or the packaging was damaged before use or the packaging is opened. Parts and tools for the implants are quite small. Pay special attention to prevent the patient from accidental aspiration of any parts during implantation!

Please do consider all recommendations in Product Catalogues issued by Callus Implant Solutions Ltd. and recommendations found in guidelines and publications of Callus Implant Systems.

Information about allergens:

Dental Implant is not recommended for patients who are allergic or sensitive for the materials listed below as the components of Callus implants:

Grade 4 Titanium: Titanium, iron (max. 0,5%), oxygen (max. 0,4%), carbon (max. 0,08%), nitrogen(max. 0,05%), hydrogen(max. 0,0125%).

Implant treatment could be considered in case of any metal sensitivity, but titanium allergy should be excluded preliminary.

MRI safety information:

Callus Implant Systems have not been evaluated for safety and compatibility in the MR environment and have not been tested for heating or migration in the MR environment. Implants are not tested in magnetic field for displacement damage, burn injury because of induction heat up or the implants shadowing effect for imaging.

Please do inform your patient after the dental implantation that he/she should warn the MR expert preliminary the MR procedure! MR diagnostic procedure could damage the implanted patient.

9. Exchange

In case the original packaging is damaged, Callus System Ltd. cannot guarantee a return implant.

10. Packaging units

In case there are no other specifications involved, implants are packed individually.

Failure to take any of the instructions detailed above excludes any and all liability of the company.

11. Signs and explanations

The signs for devices fitted with color and numerical codes can be found on the plastic handle.

Color codes

You can see it in the surgical manual.

In order to identify the correct tools and sizes please use Callus Implant System Product Catalogue!

12. Documentations and guarantee


Callus Implant Solutions Ltd recommends its users to preserve the entire documentation of each cases. The inner packaging of products contains the product code and Lot number. Inner packaging also contains a label with the product code (reference number)(REF), individual batch number (LOT) and a bar code of the product which is recommended to be documented along with patient medical data, so that each implant can be identified when necessary. The guarantee of any product is valid only with the invoice, LOT number and bar codes.



For further technical details or information please contact with your local representative!

For further information please read Callus product catalogue!

We hereby ask you if any serious incident occurs in relation to the device, that should be reported to the manufacturer and the relevant authorities of the EU Member State in which the user and/or the patient is established.

13. Explanation of the symbols on labels and instruction leaflets

 Manufacturer	 Sterilized using irradiation
 Use-by date	 Do not re-sterilize
 Date of manufacture (year, month)	 Do not use if package is damaged

 Size: Diameter (mm)	 Keep away from sunlight
 Size: length (mm)	 Keep dry
 Batch code (lot number)	 Consult instructions for use
 Catalogue number(reference number)	 Caution
 CE mark	 Rx only
 Do not reuse	United States Federal law restricts medical devices to sale by or on the order of a licensed healthcare practitioner

Distributor:
Callus Implant Solutions Ltd.
www.callusimplants.com

Manufactured:
for Callus Implant Solutions Ltd.
by Denti System Ltd.
6600 SZENTES, BESE LÁSZLÓ U. 8.
HUNGARY
version and review: 1-03/2019

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für Callus Implant Solutions Ltd.
von Denti System Ltd.